

Load Cells

By HPNE

Effortless
transportation
and weighing



Our Load Cells from AllpaQ are ideal for the easy transportation and weighing of bioprocess containers.

Created for use in a cleanroom environment, our Load Cells are ideal for the weighing of liquids during filling and discharging. These load cells can be customized for your bioprocessing needs, and the display monitor can be mounted on the unit or display mounted on a nearby wall.

Features and Benefits

- Display monitor is:
 - Mounted on the unit with a 360 deg swivel
 - Remote mounted on an adjacent wall
- Implement modifications to fit all standard and bespoke ALLpaQ units
- Connection for remote mounted display can be either via an attached cable or bluetooth
- Allow for the easy transportation and weighing of bioprocess containers in confined spaces
- Ability to be customized for logo or room label



hp-ne.com

info@hp-ne.com

401.349.4477

High Purity New England Inc.
2 Thurber Boulevard
Smithfield, RI 02917 USA